

Write a regex to create a tag group



☐ Split on underscores

☐ Data download links

Tooltip sorting method: default

Smoothing



0.6

Horizontal Axis

STEP

RELATIVE

WALL

Runs

Write a regex to filter runs

☐ ☐ test

☐ ☐ train


TOGGLE ALL RUNS

/tmp/tensorflow/mnist/logs/mnist_with_summaries

accuracy_1


1

accuracy_1



Name	Smoothed	Value	Step	Time
test	0.9336	0.9350	240.0	Tue May 2, 13:12
train	0.9284	0.9100	239.0	Tue May 2, 13:12

cross_entropy_1



dropout

1

dropout/dropout_keep_probability

